# 29

Semana  
de la **Salud**  
Ocupacional

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## SILICA, CRYSTALLINE, by XRD (filter redeposition) 7500

SiO<sub>2</sub> MW: 60.08 CAS: 14808-60-7 (quartz) RTECS: VV7330000 (quartz)  
 14464-46-1 (cristobalite) VV7325000 (cristobalite)  
 15468-32-3 (tridymite) VV7335000 (tridymite)


**METHOD:** 7500, Issue 4 **EVALUATION:** FULL **Issue 1:** 15 August 1990  
**Issue 4:** 15 March 2003

**OSHA :** quartz (respirable) 10 mg/m<sup>3</sup>/(%SiO<sub>2</sub>+2);  
 cristobalite and tridymite (respirable) ½ the above  
**NIOSH:** 0.05 mg/m<sup>3</sup>; carcinogen  
**ACGIH:** quartz (respirable) 0.1 mg/m<sup>3</sup>  
 cristobalite (respirable) 0.05 mg/m<sup>3</sup>  
 tridymite (respirable) 0.05 mg/m<sup>3</sup>


**PROPERTIES:** solid; d 2.65 g/cm<sup>3</sup> @ 0 °C; crystalline transformations: quartz to tridymite @ 867 °C; tridymite to cristobalite @ 1470 °C; α-quartz to β-quartz @ 573 °C

**SYNONYMS:** free crystalline silica; silicon dioxide

SAMPLING		MEASUREMENT	
<b>SAMPLER:</b>	CYCLONE + FILTER (10-mm nylon cyclone, Higgins-Dewell (HD) cyclone, or aluminum cyclone + 5-µm PVC membrane) *see sampling section	<b>TECHNIQUE:</b>	X-RAY POWDER DIFFRACTION
<b>FLOW RATE:</b>	Nylon cyclone: 1.7 L/min; HD cyclone: 2.2 L/min; aluminum cyclone: 2.5 L/min	<b>ANALYTE:</b>	Crystalline SiO <sub>2</sub>
<b>VOL-MIN:</b>	400 L	<b>ASH:</b>	Muffle furnace or RF plasma asher or dissolve in tetrahydrofuran
<b>-MAX:</b>	1000 L	<b>REDEPOSIT:</b>	On 0.45-µm Ag membrane filter
<b>SHIPMENT:</b>	Routine	<b>XRD:</b>	Cu target X-ray tube, graphite monochromator Optimize for intensity; 1° slit Slow step scan, 0.02°/10 sec Integrated intensity with background subtraction

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
**NIOSH:** 0.05 mg/m<sup>3</sup>; carcinogens

**OTHER OELS:** [1,2]

**PROPERTIES:** solid; d 2.65 g/cm<sup>3</sup>; crystalline transformation:  
α-quartz to β-quartz @ 573 °C; quartz to tridymite @ 867 °C;  
tridymite to cristobalite @ 1470 °C


**SYNONYMS:** free crystalline silica; silicon dioxide

SAMPLING		MEASUREMENT	
<b>SAMPLER:</b>	CYCLONE + PRE-WEIGHED FILTER (PVC), with conductive filter holder	<b>TECHNIQUE:</b>	INFRARED SPECTROMETRY
<b>FLOW RATE:</b>	~2 to ~4 L/min, specific for the cyclone used	<b>ANALYTE:</b>	Quartz / Cristobalite
<b>VOL-MIN:</b>	400 L	<b>ASHING:</b>	Muffle furnace or RF plasma asher
<b>-MAX:</b>	1000 L (total dust < 4 mg)	<b>PELLET:</b>	Mix residue with KBr; press 13-mm pellet
<b>SHIPMENT:</b>	Routine	<b>IR DETECTION:</b>	Scan absorbance between 400 cm <sup>-1</sup> and 1000 cm <sup>-1</sup>
<b>SAMPLE STABILITY:</b>	Stable	<b>CALIBRATION:</b>	Certified reference material (CRM) respirable crystalline silica (RCS) quartz & cristobalite; e.g., NIST SRM 1878b quartz, NIST 1879a cristobalite, and/or other RCS CRMs; diluted in KBr and pressed into 13-mm pellet
<b>BLANKS:</b>	3 field blanks minimum per sample set		
<b>ACCURACY</b>			

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
## FLIPPIN' AWESOME!

**ASBESTOS ALL FORMS: 0.1 f/cc (F)**

(F) Respirable fibers: length  $>5 \mu\text{m}$ ; aspect ratio  $> 3:1$ , as determined by the membrane filter method at 400-450 $\times$  magnification (4-mm objective), using phase-contrast illumination

**FUENTE: TLV'S ACGIH 2023**

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## NIOSH 7400 / OSHA ID-160

ASBESTOS and OTHER FIBERS by PCM		7400
FORMULA: Various	MW: Various	CAS: see Synonyms RTECS: Various
METHOD: 7400, Issue 2	EVALUATION: FULL	Issue 1: Rev. 3 on 15 May 1989 Issue 2: 15 August 1994
OSHA: 0.1 asbestos fiber (> 5 µm long)/cc; 1 f/cc; 30 min excursion; carcinogen	MSHA: 2 asbestos fibers/cc	NIOSH: 0.1 f/cc; fibers > 5 µm long; 400 L; carcinogen
ACGIH: 0.2 f/cc crocidolite; 0.5 f/cc amosite; 2 f/cc chrysotile and other asbestos; carcinogen	PROPERTIES: solid, fibrous, crystalline, anisotropic	
SYNONYMS (CAS #): actinolite [15669-07-5]; amosite [12172-73-5]; anthophyllite [77536-67-5]; chrysotile [12001-29-5]; serpentine [18786-24-8]; crocidolite [12001-28-4]; tremolite [77536-68-6]; amphibole asbestos [1332-21-4]; refractory ceramic fibers [143844-00-0]; fibrous glass		
SAMPLING		MEASUREMENT
<b>SAMPLER:</b> FILTER (0.45- to 1.2-µm cellulose ester membrane, 25-mm; conductive cover on cassette)	<b>TECHNIQUE:</b> LIGHT MICROSCOPY, PHASE CONTRAST	<b>ANALYTE:</b> fibers (manual count)
<b>FLOW RATE*:</b> 0.5 to 16 L/min	<b>SAMPLE PREPARATION:</b> acetone - collapse/triaceton - immersion method [2]	<b>COUNTING RULES:</b> described in previous version of this method as "A" rules [1,3]
<b>VOL-MIN*:</b> 400 L @ 0.1 fiber/cc <b>MAX*:</b> (step 4, sampling) *Adjust to give 100 to 1300 fibers/mm <sup>2</sup>	<b>EQUIPMENT:</b> 1. positive phase-contrast microscope 2. Walton-Beckett graticule (100-µm field of view) Type G-22 3. phase-shift test slide (HSE/NPL)	<b>SHIPMENT:</b> routine (pack to reduce shock)
<b>SHIPMENT:</b> routine (pack to reduce shock)	<b>CALIBRATION:</b> HSE/NPL test slide	<b>STABILITY:</b> stable
<b>BLANKS:</b> 2 to 10 field blanks per set	<b>RANGE:</b> 100 to 1300 fibers/mm <sup>2</sup> filter area	<b>BLANKS:</b> 2 to 10 field blanks per set
<b>ACCURACY</b>	<b>ESTIMATED LOD:</b> 7 fibers/mm <sup>2</sup> filter area	<b>ACCURACY</b>
<b>RANGE STUDIED:</b> 80 to 100 fibers counted	<b>PRECISION (S<sub>d</sub>):</b> 0.18 to 0.32 [1]; see EVALUATION OF METHOD	<b>RANGE STUDIED:</b> 80 to 100 fibers counted
<b>BIAS:</b> see EVALUATION OF METHOD		<b>BIAS:</b> not determined
<b>OVERALL PRECISION (S<sub>d</sub>):</b> 0.15 to 0.13 [1]		<b>OVERALL PRECISION (S<sub>d</sub>):</b> see EVALUATION OF METHOD
<b>ACCURACY:</b> see EVALUATION OF METHOD		<b>ACCURACY:</b> not determined
<b>APPLICABILITY:</b> The quantitative working range is 0.04 to 0.5 fiber/cc for a 1000-L air sample. The LOD depends on sample volume and quantity of interfering dust, and is <0.01 fiber/cc for atmospheres free of interferences. The method gives an index of airborne fibers. It is primarily used for estimating asbestos concentrations, though PCM does not differentiate between asbestos and other fibers. Fibers < ca. 0.25 µm diameter will not be detected by this method [4]. This method may be used for other materials such as fibrous glass by using alternate counting rules (see Appendix C).		
<b>INTERFERENCES:</b> If the method is used to detect a specific type of fiber, any other airborne fiber may interfere since all particles meeting the counting criteria are counted. Chain-like particles may appear fibrous. High levels of non-fibrous dust particles may obscure fibers in the field of view and increase the detection limit.		
<b>OTHER METHODS:</b> This revision replaces Method 7400, Revision #3 (dated 5/15/89).		
NIOSH Manual of Analytical Methods (NMAM), Fourth Edition		

## NIOSH 7402

ASBESTOS by TEM		7402
FORMULA: Various	MW: Various	CAS: Various RTECS: Various
METHOD: 7402	EVALUATION: PARTIAL	Issue 1: 15 May 1989 Issue 2: 15 August 1994
OSHA: 0.1 asbestos fibers (>5 µm long)/cc; 1 f/cc/30 min excursion; carcinogen	MSHA: 2 asbestos fibers/cc	NIOSH: 0.1 f/cc; fibers > 5 µm long/400 L; carcinogen
ACGIH: 0.2 crocidolite; 0.5 amosite; 2 chrysotile and other asbestos; fibrous; carcinogen	PROPERTIES: solid, fibrous, crystalline, anisotropic	
SYNONYMS (CAS #): actinolite [15669-07-5]; amosite [12172-73-5]; anthophyllite [77536-67-5]; chrysotile [12001-29-5]; serpentine [18786-24-8]; crocidolite [12001-28-4]; tremolite [77536-68-6]; amphibole asbestos [1332-21-4]		
SAMPLING		MEASUREMENT
<b>SAMPLER:</b> FILTER (0.45- to 1.2-µm cellulose ester membrane, 25-mm diameter; conductive cassette)	<b>TECHNIQUE:</b> MICROSCOPY, TRANSMISSION ELECTRON (TEM)	<b>ANALYTE:</b> asbestos fibers
<b>FLOW RATE:</b> 0.5 to 16 L/min	<b>SAMPLE PREPARATION:</b> modified Jaffe wick	<b>EQUIPMENT:</b> transmission electron microscope; energy dispersive X-ray system (EDX) analyzer
<b>VOL-MIN*:</b> 400 L @ 0.1 fiber/cc <b>MAX*:</b> (step 4, sampling) *Adjust for 100 to 1300 fibers/mm <sup>2</sup>	<b>CALIBRATION:</b> qualitative electron diffraction; calibration of TEM magnification and EDX system	<b>RANGE:</b> 100 to 1300 fibers/mm <sup>2</sup> filter area [1]
<b>SHIPMENT:</b> routine (pack to reduce shock)	<b>ESTIMATED LOD:</b> 1 confirmed asbestos fiber above 95% of expected mean blank value	<b>PRECISION (S<sub>d</sub>):</b> 0.28 when 65% of fibers are asbestos; 0.20 when adjusted fiber count is applied to PCM count [2]
<b>SHIPMENT:</b> routine (pack to reduce shock)		
<b>STABILITY:</b> stable		
<b>BLANKS:</b> 2 to 10 field blanks per set		
<b>ACCURACY</b>		
<b>RANGE STUDIED:</b> 80 to 100 fibers counted		
<b>BIAS:</b> not determined		
<b>OVERALL PRECISION (S<sub>d</sub>):</b> see EVALUATION OF METHOD		
<b>ACCURACY:</b> not determined		
<b>APPLICABILITY:</b> The quantitative working range is 0.04 to 0.5 fiber/cc for a 1000-L air sample. The LOD depends on sample volume and quantity of interfering dust, and is <0.01 fiber/cc for atmospheres free of interferences. This method is used to determine asbestos fibers in the optically visible range and is intended to complement the results obtained by phase contrast microscopy (Method 7400).		
<b>INTERFERENCES:</b> Other amphibole particles that have aspect ratios greater than 3:1 and elemental compositions similar to the asbestos minerals may interfere in the TEM analysis. Some non-amphibole minerals may give electron diffraction patterns similar to amphiboles. High concentrations of background dust interfere with fiber identification. Some non-asbestos amphibole minerals may give electron diffraction patterns similar to asbestos amphiboles.		
<b>OTHER METHODS:</b> This method is designed for use with Method 7400 (phase contrast microscopy).		
NIOSH Manual of Analytical Methods (NMAM), Fourth Edition, 8/15/94		

## OSHA ID-191

POLARIZED LIGHT MICROSCOPY OF ASBESTOS	
<b>Method number:</b>	ID-191
<b>Matrix:</b>	Bulk
<b>OSHA Content Limit:</b>	1.0%
<b>Collection Procedure:</b>	Collect approximately 1 to 2 grams of each type of material and place into separate 20 mL scintillation vials.
<b>Analytical Procedure:</b>	A portion of each separate phase is analyzed by gross examination, phase-polar examination, and central spot dispersion microscopy.
<b>Detection Limit:</b>	Less than 1% by area.
<b>Special Requirements:</b>	Send bulk samples to the laboratory in separate packages from air samples.
October 1992 Revised February 1995	Daniel T. Crane
Branch of Physical Measurements and Analysis OSHA Salt Lake Technical Center Salt Lake City, Utah 84070	
1 of 24	T-ID191-FI-01-9502-PM

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## Area Air Sampling for Asbestos Exposure

### Area Sampling Reference Guide

Sample Method*	Method Number	Flow Rate	Pump & Calibrator	Sample Media & Accessories
Asbestos by PCM (USA)	NIOSH 7400	0.5 to 16 LPM	Aircon-2 (801012-100) with power supply (801000-2 for US plug or 801000-3 for EU plug).	25mm 0.45 to 1.2 micron MCE PCM filter cassette GCU25080 (plain) or GCU25080BLG (gridded); or TEM Cassette GCS25455-50 depending upon the protocol with 1/4 inch ID Tubing
Asbestos by PCM (USA)	OSHA ID-160	0.5 to 5 LPM	Gilibrator-2 with High Flow Cell (up to 30 LPM)	25mm, 0.4 to 1.2 micron MCE
Asbestos by TEM (USA)	NIOSH 7402	0.5 to 16 LPM	GoCal (up to 20 LPM)	25mm 0.45 to 1.2 micron MCE
Asbestos by PCM (UK)	HSG248 (Replaces MDHS 39/4)	1 to 4 LPM	An optional battery pack is available (801001).	25mm diameter MCE, gridded Source UK filter cassettes locally
Asbestos by PCM (Germany) (83/477/EWG)	BGI 505-31	2 to 10 LPM	GoCal (up to 20 LPM)	37mm 0.8micron MCE Source Sampling Media through local national distributor and use FAP Holder (DEHA)
Asbestos by REM / EDXA (Germany)	BGI 505-46	2 LPM		25mm 0.4micron Polycarbonate/gold plated filter Source Sampling Media through local national distributor and use FAP Holder (DEHA)

**Background**

Asbestos is a natural mineral mined throughout the 20th century for use in fire retardant materials and vehicle brake linings. It was commonly used in buildings to improve fire safety until the 1970's, when airborne exposure to asbestos fibers was linked to lung diseases including asbestosis and lung cancer. During the mid 1980's The Occupational Health and Safety Administration in the US passed an asbestos standard that addressed worker safety in removing asbestos from buildings. Rules included air sampling for personal exposure to asbestos fibers and background & clearance sampling using area samplers.

A background sample is taken prior to the removal process, and a clearance sample is taken after the removal. During the removal process, area samples are taken inside and outside of the containment area in order to indicate the integrity of the containment area. These area samples require high volumes of air in order to attain good sensitivity. Human internal airways are covered with tiny hair-like cilia that serve to remove particulates with a wave-like motion that carries the particles from the lungs to the base of the throat, where they trigger a cough response to expel them from the body. Asbestos fibers are long and thin like microscopic spears that stick to the lung walls and cannot be removed. A secondary defense system produces tiny fibers to hold the asbestos fiber in place producing a small scar area in the lung wall representing an area that is no longer functioning. As asbestos fibers collect, the body's fibrous defense system grows rapidly and the over-production of these defense fibers is known as fibrosis; when caused by asbestos it is called asbestosis. This slow process can take decades to show up in workers. People with asbestosis are at an increased risk for lung cancer.

**Equipment for Area Sampling**

Some typical sampling methods or protocols for area asbestos sampling are listed in the table above. Some of these methods are written for both personal and area sampling and describe a wide range of flow rates. Area sampling is generally done at higher flow rates so as to produce a prescribed minimum volume of air in a reasonably short sampling time.

Analysis of the sampling media for asbestos in the various sampling methods is done by fiber count, either optically by phase contrast microscope (PCM) or by electron microscope. Various types of electron microscopes are used including transmission electron microscope (TEM), reflection electron microscope (REM) or energy dispersive X-Ray (EDXA). The sampling procedures for both optical and electron microscope laboratory methods are similar, however the filter membrane material may vary to accommodate the different fiber counting technologies.

To best select the media and sampling equipment one might first consult with the intended laboratory to determine what analysis method will be used at that lab.

\* Definitions:  
PCM: Phase Contrast Microscope  
TEM: Transmission Electron Microscope  
REM: Reflection Electron Microscope  
EDXA: Energy Dispersive X-Ray (Scanning) Electron Microscope  
MCE: Mixed cellulose ester filter

Sensidyne, LP | 1000 132th Circle North, Suite 100 | St. Petersburg, Florida 33716  
Tel: 800-451-9444 / +1 727-530-3602 | E-mail: info@sensidyne.com | Web: www.sensidyne.com

A company of the **BECKHOFF** International Group

Sample Method*	Method Number	Flow Rate	Pump & Calibrator	Sample Media & Accessories
Asbestos by PCM (USA)	NIOSH 7400	0.5 to 16 LPM	Aircon-2 (801012-100) with power supply (801000-2 for US plug or 801000-3 for EU plug).	25mm 0.45 to 1.2 micron MCE PCM filter cassette GCU25080 (plain) or GCU25080BLG (gridded); or TEM Cassette GCS25455-50 depending upon the protocol with 1/4 inch ID Tubing
Asbestos by PCM (USA)	OSHA ID-160	0.5 to 5 LPM	Gilibrator-2 with High Flow Cell (up to 30 LPM)	25mm, 0.4 to 1.2 micron MCE
Asbestos by TEM (USA)	NIOSH 7402	0.5 to 16 LPM	GoCal (up to 20 LPM)	25mm 0.45 to 1.2 micron MCE
Asbestos by PCM (UK)	HSG248 (Replaces MDHS 39/4)	1 to 4 LPM	An optional battery pack is available (801001).	25mm diameter MCE, gridded Source UK filter cassettes locally
Asbestos by PCM (Germany) (83/477/EWG)	BGI 505-31	2 to 10 LPM	GoCal (up to 20 LPM)	37mm 0.8micron MCE Source Sampling Media through local national distributor and use FAP Holder (DEHA)
Asbestos by REM / EDXA (Germany)	BGI 505-46	2 LPM		25mm 0.4micron Polycarbonate/gold plated filter Source Sampling Media through local national distributor and use FAP Holder (DEHA)

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

**David S. Rodríguez Marín, PCHI**

Business Development Manager, México

Intecon, Inc.

6590 W Rogers Circle, Suites 11 & 12

Boca Raton, FL. 33487

(561) 912 9809

[droduiguez@inteconinc.com](mailto:droduiguez@inteconinc.com)

[droduiguez@innovare-ehs.com](mailto:droduiguez@innovare-ehs.com)

**Gustavo A. Palacio Z. Qco.**

LATAM Business Development Technical  
Manager

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